

PHARMACEUTICAL ENVIRONMENTAL MONITORING 3P CULTURE MEDIA WHEN YOUR PRODUCT QUALITY MATTERS MOST



Marion LOUIS, bioMérieux

PERFORMANCES STUDY - TSA 3P, TSA 3P WITH NEUTRALISERS, COUNT-TACT 3P

DESCRIPTION

Accurate Environmental Monitoring data is crucial when making product release decisions. Culture media must be robust enough in term of development and performances to withstand the rigorous conditions in today's pharmaceutical cleanrooms.



3P™ (Pharmaceutical Proven Performances) culture media in 90mm and 55mm Count-Tact formats are recommended for the environmental monitoring within classified environments of pharmaceutical industries.


They have been developed to satisfy specific requirements encountered in pharmaceutical clean environment and were validated for high performances especially on environmental strains.

This study presents the performances of the 3P culture media allowing them to be the most comprehensive solution for environmental monitoring in classified area

GROWTH PROMOTION

The growth promotion evaluation was performed on 72 microorganisms for TSA 3P and 127 microorganisms for TSA 3P with Neutralizers and Count-Tact 3P, including pharmaceutical strains sourced from ATCC and from clean room wild isolates.

The Recovery Rates were calculated from a standard 'Gold standard' TSA not irradiated. Calibrated inoculum between 10 and 100 CFU was prepared for each organism. The plates were incubated at 20-25°C or 30-35°C according to the growth requirements of the microorganisms for an appropriate culture duration.



Excellent growth promotion characteristics confirmed by testing a wide library of pharmacopoeia and wild-type strains from clean room environments to enhance detection of any contamination during environmental monitoring.

LITERATURE REFERENCES

- ISO 14698-1 - Cleanrooms and associated controlled environments. Biocontamination control. Part 1: General principles and methods.
- ISO 18593 - Microbiology of food and animal feeding stuffs - Horizontal methods for sampling techniques from surfaces using contact plates and swabs.
- USP chapter 1116: Microbiological evaluation of cleanrooms and other controlled environments.
- EU Guidelines to Good Manufacturing Practice - Medicinal Products for Human and Veterinary Use - Annex I: Manufacture of Sterile Medicinal Products.
- Guidance for Industry Sterile Drug Products Produced by Aseptic Processing - Current Good Manufacturing Practice.


NEUTRALIZATION

TSA3P with Neutralizers and Count-Tact 3P contain neutralizing agents to inactivate any residual disinfectants present on the air, operator glove or surface to be tested.

3P performances were challenged against exposure to a range of typical and representative disinfectants. The recovery was compared to a non-irradiated TSA not exposed to disinfectants. The microorganisms tested are: *A. brasiliensis*, *C. albicans*, *P. commune*, *P. aeruginosa*, *E. coli*, *S. aureus*, *K. rhizophila*, *S. epidermidis*, *B. pumilus*, *B. subtilis* and *C. striatum*.

Different inhibition concentrations of disinfectants were determined on the media and according to the microorganisms' sensitivity to each product. 200 µL of each disinfectant solution were directly spread by calculating the recovery rate compared to an enumeration on 'Gold standard' TSA.

Incubation conditions: up to 48 hours at 30-35°C for bacteria and up to 7 days at 20-25°C for yeasts and molds.




TSA 3P with Neutralizers and Count-Tact 3P provide excellent neutralization for all these disinfectants

ISOLATOR COMPATIBILITY

Specific 3P media packaging has been designed to be resistant to Vaporous Hydrogen Peroxide (VHP) and Peracetic Acid (PA) gassing cycles as part of an isolator decontamination cycle


Peracetic Acid Decontamination

Recovery rates (%) observed before (w/o) and after (w/) a PA decontamination cycle for both format (90mm plates and Count-Tact plates)



Vaporous Hydrogen Peroxide Decontamination

Recovery rates (%) observed before (w/o) and after (w/) a VHP decontamination cycle for both format (90mm plates and Count-Tact plates)



The growth properties of the media without and with VHP and PA decontamination are similar. 3P culture media are validated with Vaporous Hydrogen Peroxide and Peracetic Acid gassing cycles for isolator compatibility.

SHELF LIFE VALIDATION

Validated with comprehensive thermal shocks during shelf-life testing to guarantee consistent performance from the beginning of the shelf-life to the date of expiry.

Stability data performed in real time have permitted to determine high product shelf-lives while guaranteeing excellent and consistent performances until their last day.

	Beginning of the shelf life			Last day of the shelf life		
	Count-Tact 3P	TSA 3P	TSA 3P w/ neutralizers	Count-Tact 3P	TSA 3P	TSA 3P w/ neutralizers
Appearance	Conform	Conform	Conform	Conform	Conform	Conform
pH	Conform	Conform	Conform	Conform	Conform	Conform
GROWTH PROMOTION TEST						
% Recovery Rate between 50 and 200% compared to a 'Gold Standard' non irradiated TSA						
<i>Candida albicans</i> ATCC 10231						
<i>Aspergillus brasiliensis</i> ATCC 1640						
<i>Pseudomonas aeruginosa</i> ATCC 9027	Conform	Conform	Conform	Conform	Conform	Conform
<i>Escherichia coli</i> ATCC 8739						
<i>Staphylococcus aureus</i> ATCC 6538						
<i>Bacillus subtilis</i> ATCC 6633						
NEUTRALISATION ACTIVITY						
% Recovery Rate between 50 and 200% compared to a 'Gold Standard' non irradiated TSA without disinfectant						
<i>Candida albicans</i> ATCC 10231		N/A (no neutralizers media)			N/A (no neutralizers media)	
<i>Staphylococcus aureus</i> ATCC 6538	Conform		Conform	Conform		Conform
<i>Bacillus subtilis</i> ATCC 6633						

Resistance to laminar flow exposure

90 mm plate formats have a 30 mL fill volume optimized for air sampling applications to reinforce resistance to dehydration effects. The 3P media are tested for both active and passive air sampling followed by worst case incubation conditions to show the media resistance to cracking and shrinkage

For the passive air sampling, plates of 3 trials batches for the TSA3P ref. 43711 / 43169 and TSA 3P with Neutralizers ref. 43811 / 43819 were left to dry for 4h30 under a vertical laminar flow and then incubated 4 days at 20-25°C followed by 3 days at 30-35°C.

No signs of dehydration at the beginning, middle and at the end of the shelf-life.

CONCLUSION

bioMérieux offers to pharmaceutical industry the most comprehensive, robust and reliable solution for environmental monitoring in clean environments.

3P™ culture media represent the state-of-the-art culture media technology from both a formulation and packaging perspective.

In addition to these proven higher performances, bioMérieux guarantees:

- Completely controlled manufacturing process ensuring consistent reliability between all batches
- Unparalleled recovery expectation 75% during routine QC
- Maximum convenience for storage and inventory management thanks to a flexible storage between 2 - 25°C, Long Shelf-life and Large monolot capacity
- Security of your crucial Environmental Monitoring samples during critical transport and incubation steps with bioMérieux's unique LOCKSURE® closure system

List of microorganisms

MOULDS	GRAM POSITIVE COCCI	GRAM POSITIVE RODS	GRAM NEGATIVE RODS	YEAST	ANEROBIC STRAINS	
<i>P. commune</i> 0909600 <i>A. brasiliensis</i> ATCC 16404 <i>P. italicum</i> ATCC 4814 <i>T. mentagrophytes</i> ATCC 9533 <i>A. coenophiala</i> ATCC 46773 <i>P. aeruginosa</i> ATCC 27952 <i>Alteromonas</i> spp. 0411781 <i>A. fermentans</i> 0411768 <i>A. fumigatus</i> 1031 <i>E. flocosum</i> 0411793 <i>Aureobasidium</i> spp. 0411787 <i>Cladosporium</i> spp. 0411782 <i>M. hiemalis</i> 8977A <i>F. acuminatum</i> 60315 <i>R. pusillus</i> 0411761 <i>S. brevicaulis</i> 0411778 <i>A. terreus</i> 0411763	<i>S. aureus</i> ATCC 6538 <i>S. epidermidis</i> 1009600 <i>K. rhizophila</i> ATCC 9341 <i>A. viscosus</i> 0501102 <i>M. spp.</i> 0306755 <i>S. nomis</i> 1009601 <i>S. haemolyticus</i> 9710129 <i>S. aureus</i> 0404752 <i>K. kristinae</i> 0406750 <i>S. pyogenes</i> 7611008 <i>M. spp.</i> 0311700 <i>M. lituseus</i> 1011600 <i>S. pyogenes</i> 19615 <i>S. saprophyticus</i> 0405751 <i>F. faecium</i> 9710124 <i>M. hiemalis</i> 8977A <i>E. gallinarum</i> 9812094 <i>M. spizizenii</i> 3456 <i>M. listeri</i> 9408011 <i>S. intermedium</i> 8202100 <i>S. epidermidis</i> 9710128 <i>E. durans</i> 5415 <i>A. terreus</i> 0411763	<i>S. haemolyticus</i> ATCC 29970 <i>S. luteus</i> 0711215 <i>S. warneri</i> 0405752 <i>S. warneri</i> 9807010 <i>B. pumilus</i> 1009602 <i>B. subtilis</i> ATCC 6633 <i>C. striatum</i> 0503059 <i>B. cereus</i> ATCC 7054 <i>S. luteus</i> 1009603 <i>A. mesenteroides</i> ATCC 23386 <i>B. licheniformis</i> 0506752 <i>L. fusiformis</i> 0102066 <i>C. haemolyticum</i> 9632770 <i>R. nestorensis</i> 0708021 <i>B. cereus</i> 0012709 <i>B. firmus</i> 0302701 <i>B. magisterum</i> 0906753 <i>B. megaterium</i> 8811015 <i>B. mucidus</i> 0611704 <i>B. pumilus</i> 0310753	<i>S. simplex</i> 0405750 <i>B. breve</i> 7846 <i>C. luteus</i> 0011707 <i>C. braakii</i> 15580 <i>C. freundii</i> 8803016 <i>E. cloacae</i> 9703036 <i>K. oxytoca</i> 13182 <i>P. vulgaris</i> 13315 <i>P. stuartii</i> 1207 <i>S. paucimobilis</i> 8509117 <i>P. oryzastrans</i> 0012704 <i>L. adhaerens</i> 9704751 <i>B. cepacia</i> 25416 <i>B. vesicularis</i> 0306754 <i>E. pneumose</i> 13883 <i>L. adhaerens</i> 9704751 <i>P. aeruginosa</i> ATCC 29027 <i>P. putida</i> 0310750 <i>R. planticola</i> 9701066 <i>A. faecalis</i> 8750 <i>R. aquatilis</i> 2637704 <i>S. gallinarum</i> 7609059 <i>E. aerogenes</i> 13048 <i>E. coli</i> 25922 <i>S. faecalis</i> 0306752 <i>A. sobria</i> 43979	<i>S. maltophilia</i> 0410754 <i>A. baumannii</i> ATCC 19506	<i>C. albicans</i> ATCC 10231 <i>G. candidum</i> 2431 <i>R. mucilaginosa</i> 8506040 <i>C. neoformans</i> ATCC 32045 <i>S. cerevisiae</i> ATCC 9763 <i>C. glabrata</i> 9911601 <i>C. sake</i> 8504059 <i>Z. cerevisiae</i> ATCC 66350 <i>C. famata</i> 8706052 <i>P. granulorum</i> LMA 94-256	<i>B. fragilis</i> ATCC 25282 <i>B. vulgatus</i> ATCC 8482 <i>C. sporogenes</i> ATCC 11437 <i>P. acnes</i> ATCC 6919 <i>C. perfringens</i> ATCC 13124 <i>P. granulosum</i> 8508092 <i>P. acnes</i> 9910602