



# IDENTIFICATION & CONFIRMATION

## List of validated methods

Reference	Identification/ confirmation Method	Principle	ISO	AOAC RI	AOAC OMA	Health Canada
				Certificate #	Certificate #	Certificate #
references according to microorganisms	API®/ID32 range	Phenotypic	Included in ISO 7218 standard for confirmation of ISO standard reference methods*  Included in ISO 16140-2 validations for confirmation of all bioMérieux pathogenic detection methods			MFHPB-20 & 30 (for Salmonella and Listeria methods confirmation)
20100 - 20160	API® 20E	Phenotypic			OMA 2009.02	
411071	VITEK® MS	Maldi-toff	Included in ISO 16140-2 validations** for confirmation of bioMérieux alternative methods	PTM 011501		
21341/21342	VITEK® 2C GN VITEK® 2C GP	Phenotypic				MFHPB-20 & 30 (for Salmonella and Listeria methods confirmation)
21345	VITEK®2 BCL	Phenotypic		PTM 080801		
21341	VITEK®2 GN	Phenotypic			OMA 2011.17	
21342	VITEK®2 GP	Phenotypic			OMA 2012.02	

\*e.g.,  
ISO 6579-1: Salmonella  
ISO 11290-1&2: Listeria  
ISO 21872-1: Vibrio  
ISO 10272-1&2: Crono  
ISO 10273: Yersinia  
ISO 10272-1&2: Campy

\*\*For:  
VIDAS® SPT  
VIDAS® Easy Salmonella  
VIDAS® Campylobacter  
CampyFood® agar detection  
CampyFood® agar enumeration  
ESIA™ One Day